

Mechanical Ventilation

Indications

- **Prolonged positive pressure ventilation**
- **Increased work of breathing**

Goals

- **Increase efficiency of breathing**
- **Increase oxygenation**
- **Improve ventilation/perfusion relationships**
- **Decrease work of breathing**

Types of Systems

- **Negative Pressure Ventilator**
 - “Iron lung”
 - Allows long-term ventilation without artificial airway
 - Maintains normal intrathoracic hemodynamics
 - Uncomfortable, limits access to patient

Types of Systems

- **Positive Pressure Ventilator**
 - Uses pressures above atmospheric pressure to push air into lungs
 - Requires use of artificial airway
 - Types
 - Pressure cycled
 - Time cycled
 - Volume cycled

Positive Pressure Ventilators

- **Pressure Cycled**
 - Terminates inspiration at preset pressure
 - Small, portable, inexpensive
 - Ventilation volume can vary with changes in airway resistance, pulmonary compliance
 - Used for short-term support of patients with no pre-existing thoracic or pulmonary problems

Positive Pressure Ventilators

- **Volume cycled**
 - **Most widely used system**
 - **Terminates inspiration at preset volume**
 - **Delivers volume at whatever pressure is required up to specified peak pressure**
 - **May produce dangerously high intrathoracic pressures**

Positive Pressure Ventilators

- **Time cycled**
 - **Terminates inspiration at preset time**
 - **Volume determined by**
 - **Length of inspiratory time**
 - **Pressure limit set**
 - **Patient airway resistance**
 - **Patient lung compliance**
 - **Common in neonatal units**

Volume-Cycled Ventilator Modes

- **Controlled Mechanical Ventilation**
 - Patient does not participate in ventilations
 - Machine initiates inspiration, does work of breathing, controls tidal volume and rate
 - Useful in apneic or heavily sedated patients
 - Useful when inspiratory effort contraindicated (flail chest)
 - Patient must be incapable of initiating breaths
 - Rarely used

Volume-Cycled Ventilator Modes

- **Assist Mode**
 - **Allows patient to control ventilator rate within limits**
 - **Inspiration begins when ventilator senses patients inspiratory effort**

Assist Mode

- **Assist/Control (A/C)**
 - Patient triggers machine to deliver breaths but machine has preset backup rate
 - Patient initiates breath--machine delivers tidal volume
 - If patient does not breathe fast enough, machine takes over at preset rate
 - Tachypneic patients may hyperventilate dangerously

Assist Mode

- **Intermittent Mandatory Ventilation (IMV)**
 - Patient breathes on own
 - Machine delivers breaths at preset intervals
 - Patient determines tidal volume of spontaneous breaths
 - Used to “wean” patients from ventilators
 - Patients with weak respiratory muscles may tire from breathing against machine’s resistance

Assist Mode

- **Synchronized Intermittent Mandatory Ventilation (SIMV)**
 - **Similar to IMV**
 - **Machine timed to delay ventilations until end of spontaneous patient breaths**
 - **Avoids over-distension of lungs**
 - **Decreases barotrauma risk**

Positive End Expiratory Pressure (PEEP)

- **Positive pressure in airway throughout expiration**
- **Holds alveoli open**
- **Improves ventilation/perfusion match**
- **Decreases FiO_2 needed to correct hypoxemia**
- **Useful in maintaining pulmonary function in non-cardiogenic pulmonary edema, especially ARDS**

Positive End Expiratory Pressure (PEEP)

DANGERS

- **High intrathoracic pressures can cause decreased venous return and decreased cardiac output**
- **May produce pulmonary barotrauma**
- **May worsen air-trapping in obstructive pulmonary disease**

Continuous Positive Airway Pressure (CPAP)

- **PEEP without preset ventilator rate or volume**
- **Physiologically similar to PEEP**
- **May be applied with or without use of a ventilator or artificial airway**
 - **Requires patient to be breathing spontaneously**
 - **Does not require a ventilator but can be performed with some ventilators**

High Frequency Ventilation (HFV)

- **Small volumes, high rates**
- **Allows gas exchange at low peak pressures**
- **Mechanism not completely understood**
- **Systems**
 - **High frequency positive pressure ventilation--60-120 breaths/min**
 - **High frequency jet ventilation--up to 400 breaths/min**
 - **High frequency oscillation--up to 3000 breaths/min**

High Frequency Ventilation (HFV)

- **Useful in managing:**
 - **Tracheobronchial or bronchopleural fistulas**
 - **Severe obstructive airway disease**
 - **Patients who develop barotrauma or decreased cardiac output with more conventional methods**
 - **Patients with head trauma who develop increased ICP with conventional methods**
 - **Patients under general anesthesia in whom ventilator movement would be undesirable**

Ventilator Settings

- Tidal volume--10 to 15ml/kg (std = 12 ml/kg)
- Respiratory rate--initially 10 to 16/minute
- FiO_2 --0.21 to 1.0 depending on disease process
 - 100% causes oxygen toxicity and atelectasis in less than 24 hours
 - 40% is safe indefinitely
 - PEEP can be added to stay below 40%
 - Goal is to achieve a $\text{PaO}_2 > 60$
- I:E Ratio--1:2 is good starting point
 - Obstructive disease requires longer expirations
 - Restrictive disease requires longer inspirations

Ventilator Settings

- **Ancillary adjustments**
 - **Inspiratory flow time**
 - **Temperature adjustments**
 - **Humidity**
 - **Trigger sensitivity**
 - **Peak airway pressure limits**
 - **Sighs**

Ventilator Complications

- **Mechanical malfunction**
 - **Keep all alarms activated at all times**
 - **BVM must always be available**
 - **If malfunction occurs, disconnect ventilator and ventilate manually**

Ventilator Complications

- **Airway malfunction**
 - Suction patient as needed
 - Keep condensation build-up out of connecting tubes
 - Auscultate chest frequently
 - End tidal CO₂ monitoring
 - Maintain desired end-tidal CO₂
 - Assess tube placement

Ventilator Complications

- **Pulmonary barotrauma**
 - **Avoid high-pressure settings for high-risk patients (COPD)**
 - **Monitor for pneumothorax**
 - **Anticipate need to decompress tension pneumothorax**

Ventilator Complications

- **Hemodynamic alterations**
 - **Decreased cardiac output, decreased venous return**
 - **Observe for:**
 - **Decreased BP**
 - **Restlessness, decreased LOC**
 - **Decreased urine output**
 - **Decreased peripheral pulses**
 - **Slow capillary refill**
 - **Pallor**
 - **Increasing Tachycardia**

Ventilator Complications

- **Renal malfunction**
- **Gastric hemorrhage**
- **Pulmonary atelectasis**
- **Infection**
- **Oxygen toxicity**
- **Loss of respiratory muscle tone**

Quick Guide to Setup

- Self check and/or Calibration as needed
- Check circuit and connections
- Set Mode: Usually “Assist/Control”
- Adjust “I” time: Usually 1 second
- Set tidal volume: 10-12 ml/kg is standard
 - May need to set “Flow” based on “I” time
- Set ventilatory rate: Adult 12-16/min

Quick Guide to Setup

- Set PEEP: std 5 cm H₂O; max 20 cm H₂O
 - Caution at 10 cm H₂O and greater
- Set “Assist/SIMV Sensitivity”: -2 cm H₂O
- Set pressure alarms
- Assess patient to confirm ventilation function
 - Monitor vital signs
 - Pulse oximetry (waveform)
 - Capnography (waveform)